

Certificate of Analysis

Product Name : LB Broth , Lennox

Cat No : L2108

Lot No : 04022127

Exp. Sep. 2024

TEST	Specifications	Results
Appearance	Off-White to light yellow powder	PASS
Solubility	20.00 g/L in \leq 30 minutes	PASS
Moisture Content	Record %	4.8 %
Particle Size (100 mesh)	Record %	NEGATIVE
pH , <i>pH was taken at room temperature (~ 22°C)</i>	6.8 - 7.20	6.0
Osmolality	Record mM/Kg	6.85
Microbial Cell Culture	E.Coli, JM101	PASS
	B. Subtilis	
Supports growth of bacterial strain	S. epidermidis	

This document has been produced electronically and is valid without signature

Date of Print : 4. 13. 2021