

★ Storage

Store desiccated at or below -20°C.

★ Contents

- Product Manual
- Daidzin, >99%

ALL PRODUCTS SOLD BY GenDEPOT ARE INTENDED FOR RESEARCH USE ONLY UNLESS OTHERWISE INDICATED. THIS PRODUCT IS NOT INTENDED FOR DIAGNOSTIC OR DRUG PURPOSE

★ Shipping Condition

Ship with ice pack.

★ Description

Daidzin is a natural organic compound in the class of phytochemicals known as isoflavones. Daidzin can be found in Japanese plant kudzu (Pueraria lobata, Fabaceae) and form soybean leaves. Daidzin is the 7-O-glucodise of daidzein. Daidzin has shown the potential for the tratment of alcohol dependency (anti -disotropic) based on animal models.

★ Solubility

Soluble in DMSO.

★ M.W.

416.38

★ Molucular Formula

C₂₁H₂₀O₉

★ Related Products

Product Name	Cat No
Albumin, (IgG, Fattu Acid and Protease Free)	A0100
AffiSelect Cox IV Loading Control Antibody	A0017
AffiSelect GAPDH Loading Control Antibody, 100ug	A0039
AffiSelect Beta-Actin Loading Control Antibody, 100ug	A0042
AffiSelect a-Tubulin Loading Control Antibody, 100ug	A0050
NP-40 Lysis Buffer (2X)	N1200
Affiselect Total Protein Extraction Solution	A0710
RIPA Cell Lysis Buffer (1X) with EDTA	R4100
Xert Protease Inhibitor Cocktail Solution (100X)	P3100
Xpert Prestained Protein Marker (6.5-240 kDa)	P8502
Xpert 2 Prestained Protein Marker (10-240 kDa)	P8503
TBS Buffer, 20X, pH 7.4	T8054
West-Ez Stripping buffer	S2100
Goat anti-Rabbit IgG(H+L)-HRP, 1mg/ml	SA002
Rabbit anti-Goat IgG(H+L)-HRP, 1mg/ml	SA007
Tween 20, Molecular Biology Grade	T9100
West-Q Pico Dura ECL Solution	W3653
West-Q Femto ECL Solution	W3700
West-Ez Blocker, 5% Non-Fat Kit	W3710