

Certificate of Analysis

Product Name : DMEM, High Glucose with L-glutamine and Sodium Pyruvate

Cat No : CM002

Lot No : 07131815

TEST	Specifications	Results
Appearance	Clear, Red Solution	PASS
Sterility		
Aerobic	Negative	NEGATIVE
Anaerobic	Negative	NEGATIVE
pH <i>pH was taken at room temperature(~ 22°C)</i>	7.0 - 7.4	7.36
Osmolarity	Record	335 mM/kg
Endotoxin <i>Endotoxin testing was performed using LAL detection method.</i>	≤ 1.0 EU/ml	< 0.10 EU/ml
Mycoplasma <i>Mycoplasma test was performed using MycoAlert® Mycoplasma Detection Assay</i>	Negative	NEGATIVE
Animal Cell Culture Cell growth promotion Viability <i>Medium was supplemented with 10% FBS and tested using BHK cells grown in a monolayer in a 5% CO₂ atmosphere.</i>	3 Passages ≥ 80%	PASS > 80%

This document has been produced electronically and is valid without signature

Date of Print : 7,18, 2018