

★ Storage

Store at -20°C

✎ Contents

- Product Manual
- Oligo (dT)₁₈ Primer
5'-d(TTTTTTTTTTTTTTTTTT)-3'

ALL PRODUCTS SOLD BY GenDEPOT ARE INTENDED FOR RESEARCH USE ONLY UNLESS OTHERWISE INDICATED. THIS PRODUCT IS NOT INTENDED FOR DIAGNOSTIC OR DRUG PURPOSE

★ Shipping Condition

Ship with ice pack.

★ Description

This primer is a synthetic single-stranded 18-mer oligonucleotide with 5'- and 3'-hydroxyl ends. The primer is supplied as a ready-to-use and can be used for first strand cDNA synthesis.

★ Usage

Applications

First strand cDNA synthesis

Concentration

100uM (0.5ug/ul)

Quality Control

Functionally tested in first strand cDNA synthesis reaction using 1.1kb RNA with 3'-poly(A) tail.

★ Related Products

Product Name	Cat No
Oligo dT (18)	O1024
PureXtract RNAsol, Trizol equivalent	R6101
amfiRivert cDNA Synthesis Platinum Master Mix	R5600
RNazor, RNase Decontamination Solution	R5800
amfiRivert Platinum One cDNA Synthesis Master Mix (4x)	R6100
RNazor, RNase Decontamination Solution	R7000
Water, DEPC Treated	W0805