

Certificate of Analysis

amfiRivert cDNA Synthesis Premix

Cat No : R5600

Lot No : 09082014

Exp. Sep. 2022

TEST	Specification	Result
Appearance	Lyophilized in 0.2ml PCR tube	PASS
DNase/RNase	None	None
Functional Test		PASS
<p>- 2-Step RT-PCR Test</p> <p>This product is tested functionally using 10 ng of total lung RNA as a templates</p> <p>1 Step - cDNA Synthesis - 42°C 30 min</p> <p>2 Step— PCR (35 cycle) using 1 ul of cDNA /GAPDH primer</p> <p>10 ul of PCR product loaded on the 1 % agarose gel and stained with SafePinky DNA staining solution for 20 min.</p> <p>Target DNA detected on UV illuminator.</p>		

This document has been produced electronically and is valid without signature

Date of Print : 9. 10. 2020