

# Certificate of Analysis

## 20X SSPE Buffer

Cat No : S0510

Lot No : 04051917

TEST	Specification	Result
Appearance	Clear, colorless liquid	PASS
pH @ 25°C : 1X Solution	7.30 - 7.40	7.34
Concentration (1X)		
Sodium Chloride	3 M $\pm$ 0.25	PASS
Sodium Phosphate	200 mM $\pm$ 0.014 $\pm$ 0.10	PASS
EDTA, Disodium Salt	20 mM $\pm$ 0.25	PASS

This document has been produced electronically and is valid without signature

Date of Print : 4.30 . 2019