

Certificate of Analysis

Goat Anti-Mouse IgG(H+L) -HRP, 200ug/500ul , HSA

Cat No : SA201

Lot No : 09051811

TEST	Result
Product Type	Secondary antibody-HRP
Source	Goat serum
E/P Ratio	The enzyme/antibody protein ration , 4:1
Form	Liquid

Specification

Tested by immunoelectrophoresis, gel diffusion and ELISA techniques as applicable. This product reacts specifically with mouse IgG and may recognize other immunoglobulin types that have light chains in common. Cross-reactivity to human serum has been minimized using affinity procedures. No antibody is detectable to non-immunoglobulin serum components.

Caution : Avoid contact with eyes, skin and clothing.

Wash all areas of contact with copious amounts of water.

See Material Safety Data Sheet for additional information.

For laboratory or further manufacturing use only.

Not intended for household use.

This document has been produced electronically and is valid without signature

Date of Print : 9. 7. 2018