

Certificate of Analysis

50X TAE Buffer

Cat No : T8050

Lot No : 02151913

| TEST | Specification | Result |
|---|---------------|--------|
| Appearance | Clear liquid | PASS |
| pH @ 25°C : 50X Solution | 8.30 – 8.40 | 8.37 |
| DNA Electrophoresis Performance Test | Tested | PASS |
| 1kb DNA Ladder separated on 1% Agarose gel. | | |

This document has been produced electronically and is valid without signature

Date of Print : 2.22 . 2019