

# Certificate of Analysis

**TBST Buffer(10X), pH 7.4**

Cat No : T8056

Lot No : 11031718

Expiration data : Nov-2021

TEST	Specification	Result
Appearance	Clear, colorless liquid	PASS
pH @ 25°C : 1X Solution	7.4 $\pm$ 0.15	PASS
Concentration (1X)		
Tris	25 mM $\pm$ 0.0025	PASS
Sodium Chloride	140 mM $\pm$ 0.014	PASS
Potassium Chloride	3.0 mM $\pm$ 0.0003	PASS
TWEEN 20	0.05% $\pm$ 0.01	PASS

This document has been produced electronically and is valid without signature

Date of Print : 11.5 . 2017